


# CURED<sup>®</sup>

HEALTH IS WHAT WE REPRESENT



CURED CATALOGUE  
2025 (ENGLISH) V.8

## Company Introduction



**Company Name:** Sepid Tan Parniyan Pouyesh Co.

**Brand:** Cured orthopedic

**Company National ID:** 10700193643

**Main Activity:** Manufacturer of Orthopedic & Physiotherapy Appliance & General Medical Devices & Supplies

**Manufacture Address:** No. 668, Farhang St, Karimabad Village, Golestan State, Iran. Postal Code: 4933147878

**Office Address:** Farhang St, Karimabad Village, Golestan State, Iran. Post Office Box: 49175-715


**Chief Executive Officer & Founder:** Seyedhamed Nabavi, WhatsApp: +989120255191

**Tel:** (+98)1732604792-3, **Cell Phone:** +989120255191, +989117939713, +989385730286, +989102710310

**Websites:** [www.curedmedical.com](http://www.curedmedical.com) / [www.curedmedical.ir](http://www.curedmedical.ir)


**Email:** [info@curedmedical.com](mailto:info@curedmedical.com)

Sepid Tan Parniyan Pouyesh Company was founded in 2013 and launched its productive activity by designing and manufacturing different types of orthosis for physiotherapy and rehabilitation in the beautiful province of Golestan in the north of Iran. Our Company have two registered CURED and XCare Brands as its product's trademark. From the beginning, this company was committed to offering high-quality products, so we worked closely with a team of orthosis and prosthesis experts in our R&D department to reach this goal till now. In 2020, after the Covid-19 epidemic worldwide, with an approach to help prevent the pandemic of this disease worldwide, we saw our duty to produce a variety of specialized medical masks such as three-layer surgical, N95, N99 and Electric Fan masks. Therefore, we have developed our R&D team and added three new production lines for producing various types of specialized medical masks and machines and Nano-filters to our production lines. On the other hand, with our self-sufficiency approach in the development of the production of specialized medical masks, we also added the production line for the production of specialized mask production machines to our production lines. In addition to all these production lines, we started a new production line of Nano-filtration fabrics by Electrospinning method with the aim of self-sufficiency in the production of specialized filtration fabrics with nanotechnology for use in medical masks, HEPA and other specialized medical filters. Today, we are self-sufficient in the manufacture of machines for the production of masks, filtration fabrics and different types of medical masks, and we can have a share in the fight against epidemic diseases such as Covid-19. Our company has a manufacturing license in three specialized groups of physiotherapy and rehabilitation and general medical devices and supplies and general medical consumption (including N95, N99, 3 Layers Surgical Mask) from the Ministry of Health and Medical Education and the Food and Drug Administration of Iran (IFDA). Also, our company is one of the knowledge-based companies located in Golestan Science and Technology Park. It is worth mentioning that, Sepid Tan Parniyan Pouyesh Co. has implemented the necessary international standards in its production line by gaining the international quality management certificates ISO13485, ISO9001 and GMP. In 2014, this company became a member of the Iranian Association of Manufacturers and Exporters of Medical, Dental and Laboratory Equipment. In order to provide services to meet the needs of customers, the investigation and feasibility of providing up-to-date products, under the supervision of the experienced experts, with the aim of high-quality manufacturing have been done. Our R & D department, after research & reviewing the performance of different devices and equipment, has succeeded in manufacturing different types of medical equipment and their production line machines. Today, combining cutting-edge technology and highly professional staff also innovative ideas along with competitive prices, highly focused R&D and efficient communication, have enabled the company to maintain its top position in the Middle East.




# Spine

The spine has three main sections, cervical, thoracic, and lumbar. There are seven cervical vertebrae. The twelve thoracic and five lumbar vertebrae make up the thoracolumbar spine. The cervical spine supports the head and allows anterior, posterior, and lateral flexion and extension, plus rotation of the head and neck. C2, called the axis, forms a pivot with C1, called the atlas, around which the head rotates on the cervical spine. Normal curvature is lordotic (eyes pointing up). The thoracolumbar spine (T1-L5) has two normal curvatures, lordotic in the lumbar region and kyphotic in the thoracic region. The primary spinal functions are to maintain an erect posture, provide stability and mobility, transmit loads, absorb shocks, and protect the spinal cord.



# Spine Cervical



## Soft Collar Cervical Collar

Comfortable  
Protection &  
Covered with two  
layers of fabric  
Code: 1201

Comfortable  
Protection



Covered with two  
layers of fabric



### Indications:


- Cervical spine contusion
- Acute wry neck
- Cervical spondylosis
- Cervical strain/sprain
- HIVD (Herniated Intervertebral Disc)
- Post-operative

### Features:

- Cervical collar, anatomically shaped. Hook and loop closure allows for easy application.
- Two layers fabric covered.
- Supports the cervical spine in a neutral position.
- Available in three heights for adults version: 7.5cm (2.95 inch) 9cm (3.54 inch), 10.5cm (4.13 inch)
- Available in three heights for Kids version: 4.5cm (1.7 inch) 5.5cm (2.16 inch), 6.5cm (2.5 inch)

### Table of Size:

- To measure the height of the neck, eyes should be looking straight ahead then measure the distance between the edge of sternum and lower jaw with the neck
- To measure the circumference of the neck, measure from the circumference of the neck at "Adam's" apple.



S

32-36cm (12.59-14.17")

M

36-40cm (14.17-15.74")

L

40-44cm (15.74-17.32")

XL


44-48cm (17.32-18.89")

2XL

48-52cm (18.89-20.47")

### Materials:

TDI Foam, Elastane, Nylon, Polyester



# Spine Cervical



Comfortable Protection

Covered with two layers of fabric

Modular Design

## Semi Rigid Splint

### Cervical Collar

Comfortable Protection & Covered with two layers of fabric  
Code: 1202

#### Indications:

- Acute wry neck
- Cervical strain/sprain
- HIVD (Herniated Intervertebral Disc)
- Arthritis
- Post-operative
- Mild scoliosis


#### Features:

- Cervical collar, anatomically shaped.
- Two layers fabric covered.
- Provides adjustable immobilization of the cervical spine and supports the cervical spine in a neutral position.
- Soft foam construction with integrated reinforcement splint provides firm support for the unstable neck.
- The vinyl panel can be detached to allow for a gradual increase in range of motion.
- Velcro Strap allows for easy application.
- Available in three heights: 7.5cm (2.95 inch) 9cm (3.54 inch), 10.5cm (4.13 inch)

#### Table of Size:

- To measure the height of the neck, eyes should be looking straight ahead then measure the distance between the edge of sternum and lower jaw with the neck in the extension position.

- To measure the circumference of the neck, measure from the circumference of the neck at "Adam's" apple.




S	32-36cm (12.59-14.17")
M	36-40cm (14.17-15.74")
L	40-44cm (15.74-17.32")
XL	44-48cm (17.32-18.89")
2XL	48-52cm (18.89-20.47")

#### Materials:

TDI Foam, Polyethylene, Elastane, Nylon, Polyester

# Spine Cervical



## Hard ADJ Cervical Collar

Comfortable  
Protection  
Code: 1301



Comfortable  
Protection



### Indications:

- HIVD (Herniated Intervertebral Disc)
- Post-operative
- Cervical fracture
- Cervical spondylosis

### Features:


- Supports the cervical spine in a neutral position.
- Cervical collar, anatomically shaped, foam edge padding provides enhanced cushioning and comfort. It ensures long functional life. It is extremely light in weight.
- Anatomically designed with height adjustment ensures perfect fitting, improves immobilization and enhances comfort.
- High-quality PE Sheet ensures rigid support, long life of the collar and good aesthetics.
- Hook and loop closure allows for easy application.
- The perforated body improves ventilation, enhances comfort and compliance.
- This is a height-adjustable cervical collar, also it is available in three heights: 7.5cm (2.95 inches), 9cm(3.54 inches), 10.5cm (4.13 inches) for all different neck heights.

### Materials:

Polyethylene, Foam, Nylon, Polyester.


### Table of Size:

- To measure the height of the neck, eyes should be looking straight ahead then measure the distance between the edge of sternum and lower jaw with the neck in the extension position.
- To measure the circumference of the neck, measure from the circumference of the neck at "Adam's apple."



**Free Size** 30-44 cm (11.81-17.32")

**Over Size** 34-48 cm (13.38-18.89")



# Spine Cervical



Used in Hospitals



## Philadelphia Cervical Collar

Used in Hospitals

Code: 1101  
(Tracheostomy & MRI  
Compatible)

Code: 1102  
(MRI Compatible)

### Indications:

- Temporary fixation after cervical spine injury
- After treatment of cervical spine after injuries and operations
- Discopathy of the area of cervical spine in want of temporary fixation
- Acute torticollis

### Features:


The CURED Philadelphia Cervical Collars produce in three design that is shown in the top picture. These products are constructed from Foam and provide firm support designed to help immobilize the neck following soft tissue or disc trauma. Especially useful following surgery and in rehabilitation. Code 1101 and 1102 have MRI compatibility, Code 1101 also has an opening to allow tracheal access. The moulded structure helps place the cervical spine in slight protraction, helping with relief for flexion injuries and promoting correct cervical posture and alignment. This collar helps cervical musculature to rest and recover by reducing movement at the neck and ensuring muscles and tendons are not overly strained. It is designed to mould comfortably above the

collarbones and under the chin, helping to relieve some gravitational pressure caused by the weight of the head on the neck.

### Table of Size:

- Eyes should be looking straight ahead.
- Measure distance between the edge of sternum and lower jaw with the neck in the extension position.

XS	28-32 cm (11.02-12.59")
S	32-36 cm (12.59-14.17")
M	36-40 cm (14.17-15.74")
L	40-44 cm (15.74-17.32")
XL	44-48 cm (17.32-18.89")



### Materials:

Foam, EVA, PE, Nylon, Polyester.



# Spine Thoracic



## Posture Brace

w/Steel Stay

Two  
Flexible Stays  
General: 3101  
Economic: 31012



### Indications:

- Early kyphosis
- Poor posture
- Weak scapula muscles
- Round shoulder

### Features:


- Additional flexible stays help reinforce back alignment and proper posture and prevent slouching.
- It can also help ease muscular aches and spasms by supporting the soft tissues and evenly distributing forces over the lumbar area.
- Triangle configuration correctly distributes pressure to the tip of the shoulder.
- Crisscrossing straps provide more support and make the brace easy to adjust also provide a custom fit.
- Slimline, so it can be worn under everyday clothes.
- Kids version available.

### Materials:


Neoprene, Cotton, Nylon, Polyester, Steel.

### Table of Size:

Circumference of chest measured below the breast.




A = Men: Circumference of chest measured above breast



B = Women: Circumference of chest measured below breast

<b>S</b>	64-74 cm (25.20-29.13")	<b>XL</b>	94-104 cm (37.01-40.94")
<b>M</b>	74-84 cm (29.13-33.07")	<b>2XL</b>	104-114 cm (40.94-44.88")
<b>L</b>	84-94 cm (33.07-37.01")	<b>3XL</b>	114-124 cm (44.88-48.81")



# Spine Thoracic



## Clavicle Brace

Adjustable Fit  
Code: 3301



### Indications:

- Clavicle fracture
- Post- surgical support
- Poor posture

### Features:


- Clavicle brace is constructed from lightweight and durable thick foam and breathable material.
- Has considerable wide shoulder straps with generous padding that, provide proper immobilization and enhances patient comfort.
- Has thick pad at back helps to prevent irritation.
- Has hook and loop fasteners for easy fixing.
- The top strap can be adjusted to manage the rigidity of fixation by down to the up movement.
- Movable metal ring for proper fit.
- Suitable for adults or children.

### Materials:


Cotton, Nylon, Polyester, Foam, Metal.

### Table of Size:

Measure around the chest.



A = Men: Circumference of chest measured above breast



B = Women: Circumference of chest measured below breast

S

65-80 cm (25.59-31.49")

M

80-95 cm (31.49-37.40")

L

95-110 cm (37.4-43.30")

XL

110-125 cm (43.30-49.21")

2XL

125-140 cm (49.21-55.11")


# Spine & Abdominal



## Abdominal Binder

w/Compression Strap

Used in Hospitals  
Code: 4202



### Indications:

- After birth
- Mild low back pain
- Muscle strain/ sprain
- After abdominal operation

### Features:

- Supports lower torso and reduces muscle fatigue and strain.
- Ideal for post-surgical use for effective support without the stiffness of a corset.
- Designed to lift and trim the abdomen, the abdomen support fits discreetly under clothing.
- Breathable elastic mesh for great air-flow to reduce sweating.
- Double elastic straps with hook and loop closure for easy size adjustment, allows support adjustable from gentle, to strong compression.
- Compression strap provides extra compression.
- The height of the abdominal part is 250mm.
- Universal design is intended for both men and women.
- Free Size is available.


### Materials:

Elastane, Nylon, Polyester.

P10

### Table of Size:

Measure around the waist.



<b>S</b>	70-80 cm (27.56-31.50")	<b>2XL</b>	110-120 cm (43.31-47.24")
<b>M</b>	80-90 cm (31.50-35.43")	<b>3XL</b>	120-130 cm (47.24-51.18")
<b>L</b>	90-100 cm (35.43-39.37")	<b>4XL</b>	130-140 cm (51.18-55.12")
<b>XL</b>	100-110 cm (39.37-43.31")	<b>5XL</b>	140-150cm (55.12-59.06")
<b>F</b>	70-120 cm (27.56-47.24")		



# Spine & Abdominal



## Neoprene Waist Trimmer

Code: 4101



### Indications:

- Helps support the abdominal region.
- Muscle strain/sprain.
- Mild low back pain.
- Post-surgical abdominal support.

### Features:


- Helps support the abdominal region.
- Easily adjusts for a comfortable fit.
- Designed to lift and trim the abdomen.
- Universal design is intended for both men and women.

### Materials:

Neoprene, Nylon, Elastane, cotton, Polyester.

### Table of Size:

Measure around the waist.



S	70-80 cm (27.56-31.50")
M	80-90 cm (31.50-35.43")
L	90-100 cm (35.43-39.37")
XL	100-110 cm (39.37-43.31")
2XL	110-120 cm (43.31-47.24")
3XL	120-130 cm (47.24-51.18")
4XL	130-140 cm (51.18-55.12")
5XL	140-150cm (55.12-59.06")

# Spine Lumbar



## Low Profile Industrial Back Support

w/Shoulder  
Straps

Code: 5203



### Indications:

- Low back pain
- Muscle strain/ligament sprain
- Prevent injury

### Features:


- Helps prevent and relieve lower back pain.
- Easy to apply and remove. Inconspicuous to wear.
- May be used for both the lumbar and lumbosacral regions.
- Four pliable support splints in the back for extra stability.

### Materials:

Nylon, Polyester, Elastane, PE.


### Table of Size:

Measure around the waist.



S	70-80 cm (27.56-31.50")
M	80-90 cm (31.50-35.43")
L	90-100 cm (35.43-39.37")
XL	100-110 cm (39.37-43.31")
2XL	110-120 cm (43.31-47.24")
3XL	120-130 cm (47.24-51.18")
4XL	130-140 cm (51.18-55.12")
5XL	140-150cm (55.12-59.06")

# Spine Lumbar



## Sacro Lumbar Belt w/Power Straps

Four Splints  
Code: 5204(27cm)  
Code: 5205(30cm)

### Indications:

- Spondylosis
- Low back pain
- Spondylolisthesis
- Muscle strain/ sprain
- HIVD (Herniated Intervertebral Disc)

### Features:


- The 5204 model has a 27cm back height and the 5205 model has a 30cm back height.
- Special breathable material for a comfortable fit.
- Four posterior flexible stays stabilize the spine, allowing limited flexion and extension.
- Power straps help provide added support, compression and serve to increase rigidity and comfort to tighten the appliance if needed.
- The height of the belt in the front part is shortened.
- Four splints in back for stabilization and to help prevent rolling.

### Materials:

Elastane, Nylon, Polyester, Polyethylene.

### Table of Size:

Measure around the waist.



S	70-80 cm (27.56-31.50")
M	80-90 cm (31.50-35.43")
L	90-100 cm (35.43-39.37")
XL	100-110 cm (39.37-43.31")
2XL	110-120 cm (43.31-47.24")
3XL	120-130 cm (47.24-51.18")
4XL	130-140 cm (51.18-55.12")
5XL	140-150cm (55.12-59.06")

# Spine Lumbar



## Sacro Lumbar Belt

w/Power Straps  
(Econ.)

Four Splints  
Code: 52042(27cm)  
Code: 52052(30cm)



### Indications:

- Spondylosis
- Low back pain
- Spondylolisthesis
- Muscle strain/ sprain
- HIVD (Herniated Intervertebral Disc)

### Features:


- The 52042 model has a 27cm back height and the 52052 model has a 30cm back height.
- Special breathable material for a comfortable fit.
- Four posterior flexible stays stabilize the spine, allowing limited flexion and extension.
- Power straps help provide added support, compression and serve to increase rigidity and comfort to tighten the appliance if needed.
- The height of the belt in the front part is shortened.
- Four splints in back for stabilization and to help prevent rolling.

### Materials:

Elastane, Nylon, Polyester, Polyethylene.

### Table of Size:

Measure around the waist.



S	70-80 cm (27.56-31.50")
M	80-90 cm (31.50-35.43")
L	90-100 cm (35.43-39.37")
XL	100-110 cm (39.37-43.31")
2XL	110-120 cm (43.31-47.24")
3XL	120-130 cm (47.24-51.18")
4XL	130-140 cm (51.18-55.12")
5XL	140-150cm (55.12-59.06")

# Spine Lumbar



## Sacro Lumbar Belt

w/Power Strap  
(Free Size)

Four Splints  
Code: 52041(27cm)  
Code: 52051(30cm)

### Indications:

- Spondylosis
- Low back pain
- Spondylolisthesis
- Muscle strain/ sprain
- HIVD (Herniated Intervertebral Disc)

### Features:


- The 52041 model has a 27cm back height and the 52051 model has a 30cm back height.
- Special breathable material for a comfortable fit.
- Four posterior flexible stays stabilize the spine, allowing limited flexion and extension.
- Power straps help provide added support, compression and serve to increase rigidity and comfort to tighten the appliance if needed.
- The height of the belt in the front part is shortened.
- Four splints in back for stabilization and to help prevent rolling.
- This product is One Size.

### Materials:

Elastane, Nylon, Polyester, Polyethylene.

### Table of Size:

Measure around the waist. This product is one size that fits most waist support and will fit.



F

70-120 cm (27.56-47.24")

# Spine Lumbar



## Sacro Lumbar Belt

w/Power Straps  
(Free size & Econ.)

Four Splints  
Code: 520412(27cm)  
Code: 520512(30cm)



### Indications:

- Spondylosis
- Low back pain
- Spondylolisthesis
- Muscle strain/ sprain
- HIVD (Herniated Intervertebral Disc)

### Features:


- The 520412 model has a 27cm back height and the 520512 model has a 30cm back height.
- Special breathable material for a comfortable fit.
- Four posterior flexible stays stabilize the spine, allowing limited flexion and extension.
- Power straps help provide added support, compression and serve to increase rigidity and comfort to tighten the appliance if needed.
- The height of the belt in the front part is shortened.
- Four splints in back for stabilization and to help prevent rolling.
- This product is One Size that Fits Most waist support and will fit.

### Materials:

Elastane, Nylon, Polyester, Polyethylene.


### Table of Size:

Measure around the waist



F

70-120 cm (27.56-47.24")



# Spine Lumbar



1:6 Pulley System  
**6X**

Adjustable Fit

## General LSO Flex

w/1:6 Pulley System  
Compression

Corset offers  
6:1compression  
and control  
from L1 to S1.  
Code: 5301

### Indications:

- Spondylosis
- Low back pain
- Sprains and strains
- Spondylolisthesis

### Features:


- Advanced mechanical 6:1 pulley system provides compression, control of spine vertebrae, and intra-abdominal compression with minimal effort.
- Single pull-strap closure provides easy and quick one-handed adjustment.
- Patient-friendly, especially for arthritic patients.
- Washable and metal-free construction.
- Medical Imaging friendly.
- Free Size is available.

### Materials:

Cotton, Nylon, Polyester, Polyethylene


### Table of Size:

Measure around the waist.



<b>S</b>	70-80 cm (27.56-31.50")	<b>2XL</b>	110-120 cm (43.31-47.24")
<b>M</b>	80-90 cm (31.50-35.43")	<b>3XL</b>	120-130 cm (47.24-51.18")
<b>L</b>	90-100 cm (35.43-39.37")	<b>4XL</b>	130-140 cm (51.18-55.12")
<b>XL</b>	100-110 cm (39.37-43.31")	<b>5XL</b>	140-150cm (55.12-59.06")
<b>F</b>	70-130 cm (27.56-51.18")		

# Spine Lumbar



## Pro LSO Flex

w/1:6 Pulley System  
Comp. 23/30/35/45 cm



Corset offers 6:1  
compression and  
control from L1 to S1.  
Code: 5302



### Indications:

- Low back pain.
- Sprains and strains.
- Lumbar disc displacement, osteoporosis.
- Disc degeneration and herniation.
- Postoperative fusion, laminectomy, discectomy.
- IDET procedures.
- Spondylolisthesis.
- Spinal stenosis.

### Features:

- Advanced mechanical 6:1 pulley system provides compression, control of spine vertebrae, and intra-abdominal compression with minimal effort.
- Single pull-strap closure provides easy and quick one-handed adjustment.
- Patient-friendly, especially for arthritic patients.
- Modular design, progressive support & control.
- Washable and metal-free construction.
- Rigid posterior removable extra back panel increases support.
- The anatomically contoured posterior panel can be heat moulded and customized to fit individual patients.
- Padded panels increase patient comfort.


- Available in Four back panel heights: 23cm(9.06 inches), 30cm(11.81 inches), 35cm(13.77 inches), 45cm(17.7 inches).
- Free Size available.
- Medical Imaging friendly.

### Materials:


Cotton, Nylon, Polyester, Polyethylene.

### Table of Size:

Measure around the waist.



S	70-80 cm (27.56-31.50")	2XL	110-120 cm (43.31-47.24")
M	80-90 cm (31.50-35.43")	3XL	120-130 cm (47.24-51.18")
L	90-100 cm (35.43-39.37")	4XL	130-140 cm (51.18-55.12")
XL	100-110 cm (39.37-43.31")	5XL	140-150cm (55.12-59.06")
F	70-130 cm (27.56-51.18")		



# Spine Lumbar



## Expert LSO Flex

w/1:6 Pulley System  
Comp. 23/30/35 cm



6X



Corset offers 6:1  
compression and  
control from L1 to S1.  
Code: 5303

### Indications:

- Low back pain.
- Sprains and strains.
- Lumbar disc displacement, osteoporosis.
- Disc degeneration and herniation.
- Postoperative fusion, laminectomy, discectomy.
- Spondylolisthesis.
- Spinal stenosis.

### Features:

- Advanced mechanical 6:1 pulley system provides compression, control of spine vertebrae, and intra-abdominal compression with minimal effort.
- Single pull-strap closure provides easy and quick one-handed adjustment.
- Patient-friendly, especially for arthritic patients.
- Modular design, progressive support and control.
- Washable and metal-free construction.
- Rigid posterior removable extra back panel increases support.
- The anatomically contoured posterior panel can be heat moulded and customized to fit individual patients.
- Free Size available.


- Lateral pads increase support and prevent pressure points.
- Padded panels increase patient comfort.
- Available in three back panel heights: 23cm(9.06 inches), 30cm(11.81 inches), 35cm(13.77 inches).
- Medical Imaging friendly.

### Materials:

Cotton, Nylon, Polyester, Polyethylene.


### Table of Size:

Measure around the waist.



S	70-80 cm (27.56-31.50")	2XL	110-120 cm (43.31-47.24")
M	80-90 cm (31.50-35.43")	3XL	120-130 cm (47.24-51.18")
L	90-100 cm (35.43-39.37")	4XL	130-140 cm (51.18-55.12")
XL	100-110 cm (39.37-43.31")	5XL	140-150cm (55.12-59.06")
F	70-130 cm (27.56-51.18")		

# Spine Lumbar



## Expert LSO

w/1:6 Pulley System  
Compression 40cm

Corset offers 6:1  
compression and  
control from T9 to S1  
Code: 5304



1:6 Pulley  
System

6X

Adjustable  
Fit & Height



Prescription  
Only



### Indications:

- Lumbar back pain.
- Lumbar strains/sprains.
- Lumbar spinal stenosis.
- Disc protrusions.
- Spondylosis /Spondylolisthesis.
- Disc degeneration and herniation.
- Postoperative fusion, laminectomy, discectomy.
- Pre- & Post-surgical treatment

### Features:

- The product features a rigid posterior panel that extends from the sacrococcygeal junction up to the T9 level (spine of the scapula).
- Advanced mechanical 6:1 pulley system provides compression, control of spine vertebrae, and intra-abdominal compression with minimal effort.
- Single pull-strap closure provides easy and quick one-handed adjustment.
- Patient-friendly, especially for arthritic patients.
- Modular design, progressive support & control.
- Washable and metal-free construction.
- Rigid posterior removable extra back panel increases support.
- The anatomically contoured posterior panel can be heat moulded and customized to fit individual patients.


- Padded panels increase patient comfort.
- Lateral and abdominal pads increase support and prevent pressure points.
- Available in 40cm(15.7 inches) back panel height
- Free Size available.

### Materials:


Cotton, Nylon, Polyester, POM, PE, PP, PVC.

### Table of Size:

Measure around the waist.



S	70-80 cm (27.56-31.50")	2XL	110-120 cm (43.31-47.24")
M	80-90 cm (31.50-35.43")	3XL	120-130 cm (47.24-51.18")
L	90-100 cm (35.43-39.37")	4XL	130-140 cm (51.18-55.12")
XL	100-110 cm (39.37-43.31")	5XL	140-150cm (55.12-59.06")
F	70-130 cm (27.56-51.18")		



# Spine Thoracic & Lumbar



1:6 Pulley System

**6X**

Adjustable Fit & Height

Prescription Only

## Expert TLSO w/1:6 Pulley System Compression

Corset offers 6:1 compression and control from T3 to S1  
Code: 5405

### Indications:

- Spinal stenosis.
- Excessive kyphosis.
- Acute thoracic/lumbar strains.
- Spondylosis /Spondylolisthesis.
- Disc degeneration and herniation.
- Postoperative fusion, laminectomy, discectomy.
- Pre- & Post-surgical treatment

### Features:

- The CURED Expert TLSO w/1:6 advanced mechanical Pulley System Compression is designed to provide effective thoracic lumbar motion restriction and premium patient comfort.
- The product features a rigid posterior panel that extends from the sacrococcygeal junction up to the T3 level (spine of the scapula).
- Single pull-strap closure provides easy and quick one-handed adjustment.
- Modular design, progressive support and control and washable construction.
- The anatomically contoured posterior panel increases support and can be heat moulded and customized to fit individual patients.
- Padded panels increase patient comfort.
- Lateral and abdominal pads increase support and prevent pressure points.


Cured | [www.curedmedical.com](http://www.curedmedical.com)

- Adjusting the superior section of the Expert TLSO allows for modification of the orthosis, enabling its conversion into an LSO configuration.
- Available in a free size version that fits most.
- The height of the Expert TLSO is adjustable according to the patient's height.

### Materials:


Cotton, Nylon, Polyester, POM, PE, PP, PVC.

**Table of Size:** Measure around the waist.



<b>S</b>	70-80 cm (27.56-31.50")	<b>2XL</b>	110-120 cm (43.31-47.24")
<b>M</b>	80-90 cm (31.50-35.43")	<b>3XL</b>	120-130 cm (47.24-51.18")
<b>L</b>	90-100 cm (35.43-39.37")	<b>4XL</b>	130-140 cm (51.18-55.12")
<b>XL</b>	100-110 cm (39.37-43.31")	<b>5XL</b>	140-150 cm (55.12-59.06")
<b>F</b>	70-130 cm (27.56-51.18")		

# Spine Lumbar



## Maternity Support

Code: 6201



### Indications:

- Relieves lumbar, leg and abdomen muscles pain associated with pregnancy
- Back discomfort during pregnancy

### Features:


- Provides support to the lumbar and pelvic region.
- Relieves the strain on muscles ligaments caused by the pendulous abdomen.
- Helps ease back pain and provides abdominal support.
- Comfortable against the skin.
- Helps relieve pain by gently lifting the abdomen.
- Smooth and invisible under clothing.

### Materials:

Nylon, Polyester, Elastane.

### Table of Size:

Measure around the hip.




Free Size

60-110 cm (23.62-43.30")

Over Size

70-130 cm (27.55-51.18")



# Spine Lumbar



## Maternity Support Belt

Code: 6202

### Indications:

- Before and after pregnancy use.
- Instant relief from lower back pain.
- Relieves lumbar, leg, pelvic region and abdomen muscles pain associated with pregnancy.
- Helps to balance the body weight and the strain on muscles ligaments by gently lifting the extended abdomen and transferring some of its weight from front over the pelvis to back of the spine and secure pendulous abdomen side effects during pregnancy.

### Features:


- Helps relieve pain by gently lifting the abdomen
- Additional diagonal belt ensures that the pregnancy support belt is in proper position.
- Constructed in two parts and fully adjustable to accommodate abdominal growth throughout pregnancy
- Comes with 2 stays which provide extra support in lumbar region and hold the support belt in place and prevent it from rolling.
- Provides gentle support with maximum comfort without undesirable pressure.

Cured | [www.curedmedical.com](http://www.curedmedical.com)

- Comfortable and inconspicuous when in use.

### Table of Size:

Measure around the hip.



S	57-93 cm (22.44-36.61")
M	67-106 cm (26.37-41.73")
L	77-119 cm (30.31-46.85")
XL	87-132 cm (34.25-51.96")
2XL	97-145 cm (38.18-57.08")
F	61.5-118.5cm (24.21-46.65")

### Materials:

Nylon, Polyester, Elastane, PE.

# Upper Extremity



# Upper Extremity

The upper limbs or upper extremities are the forelimbs of an upright-postured tetrapod vertebrate, extending from the scapulae and clavicles down to and including the digits, including all the musculatures and ligaments involved with the shoulder, elbow, wrist and knuckle joints. In humans, each upper limb is divided into the arm, forearm and hand, and is primarily used for climbing, lifting and manipulating objects.





# Shoulder Support



## Arm Sling

Breathable  
Fabric  
Code: 7331

### Indications:

- Brachial nerve injury
- Dislocated shoulder
- Upper arm fracture
- Post-surgical protection

### Features:


- Built-in thumb loops support your hand to prevent fatigue.
- Supports the arm and carries the weight across the back and shoulder.
- Universal design fits right or left arm.
- Adjustable shoulder strap with comfort pad provides a customisable fit.
- Comfortable lightweight material.
- Kids version available.

### Materials:

Cotton, Polyester, Nylon, Polyethylene.

### Table of Size

Measure the length of forearm.



XS

15.2–20.3cm (6 – 8")

S

20.3–25.4cm (8–10")

M

25.4–30.5cm (10–12")

L

30.5–35.6cm (12–14")

XL

35.6–39.4cm (14–15.5")

# Shoulder Support



## Breathable Arm Sling

Breathable Fabric  
Code: 73315



### Indications:

- Brachial nerve injury
- Dislocated shoulder
- Upper arm fracture
- Post-surgical protection

### Features:


- Built-in thumb loops support your hand to prevent fatigue.
- Supports the arm and carries the weight across the back and shoulder.
- Universal design fits right or left arm.
- Adjustable shoulder strap with comfort pad provides a customisable fit.
- Comfortable lightweight material.

### Materials:

Cotton, Polyester, Nylon, Polyethylene.

### Table of Size

Measure the length of forearm.



<b>XS</b>	15.2–20.3cm (6 – 8")
<b>S</b>	20.3–25.4cm (8–10")
<b>M</b>	25.4–30.5cm (10–12")
<b>L</b>	30.5–35.6cm (12–14")
<b>XL</b>	35.6–39.4cm (14–15.5")
<b>F</b>	25.4–39.4cm (10–15.5")



## Upper Extremity



## Ultra Arm Sling

Breathable  
Fabric  
Code: 7302



### Indications:

- Brachial nerve injury
- Dislocated shoulder
- Upper arm fracture
- Rotator cuff injuries
- Post-surgical protection

### Features:


- Helps immobilize and support the arm.
- Used to comfortably support a sprained or broken arm or for post-surgical support.
- Built-in thumb loops support your hand to prevent fatigue.
- Adjustable shoulder strap with comfort pad provides a customisable fit.
- The removable extra strap offers more stabilization to the arm and body.
- Comfortable lightweight and breathable material.
- Universal design fits right or left arm.

### Materials:

Polyester, Nylon, Polyethylene.

### Table of Size:

Measure the length of forearm.



<b>XS</b>	15.2–20.3cm (6 – 8")
<b>S</b>	20.3–25.4cm (8–10")
<b>M</b>	25.4–30.5cm (10–12")
<b>L</b>	30.5–35.6cm (12–14")
<b>XL</b>	35.6–39.4cm (14–15.5")

# Shoulder Support



## Universal Shoulder Immobilizer

Adjustable Fit  
Code: 7303



### Indications

- Brachial nerve injury
- Dislocated shoulder
- Upper arm fracture
- Rotator cuff injuries
- Impingement pain
- Post-surgical protection

### Features:


- Secures arm against the body and restricts the mobility of the shoulder and upper arm for separations, dislocations, strains and other injuries.
- Adjustable hook and loop closures allow swathe to be positioned in the desired location.
- May be used as an arm sling if swathe is not desired.
- Adjustable straps provide a customisable fit.
- Comfortable lightweight material.
- Universal design fits right or left arm.

### Materials:


Cotton, Nylon, Polyester, PE.

### Table of Size:

Measure around the chest.



A = Men: Circumference of chest measured above breast



B = Women: Circumference of chest measured below breast

S

60.96-91.44cm (24-36")

M

76.20-106.68cm (30-42")

L

91.44-121.92cm (36-48")

XL


106.68-137.16cm (42-54")

F

76-137.16cm (30-54")



# Shoulder Support



## Universal Sling & Swathe

Adjustable Fit  
Code: 73031



### Indications

- Brachial nerve injury
- Dislocated shoulder
- Upper arm fracture
- Rotator cuff injuries
- Impingement pain
- Post-surgical protection

### Features:


- Secures arm against the body and restricts the mobility of the shoulder and upper arm for separations, dislocations, strains and other injuries.
- Adjustable hook and loop closures allow swathe to be positioned in the desired location.
- May be used as an arm sling if swathe is not desired.
- Adjustable straps provide a customisable fit.
- Comfortable lightweight material.
- Universal design fits right or left arm.
- One size product is available.

### Materials:

Cotton, Nylon, Polyester, PE.

### Table of Size:

Measure the length of forearm.



F	25.4-44.5cm (10-17.5")
XS	15.2-20.3cm (6 - 8")
S	20.3-25.4cm (8-10")
M	25.4-30.5cm (10-12")
L	30.5-35.6cm (12-14")
XL	35.6-39.4cm (14-15.5")

# Shoulder Support



## Shoulder Abduction Orthosis

Code: 7308



### Indications

- As an orthopedic positioning method for sprains, luxations and glenohumeral subluxations.
- Capsular displacement.
- As a conservative treatment for acute and chronic bursitis.
- In rotator cuff repair surgery.
- In shoulder endoprosthesis implantation.
- Adhesive capsulitis.
- Entrapment syndrome.
- As a method of immobilization after implantation of a prosthetic the shoulder joint and after acromioplasty.


### Features:

- Manufactured in breathable honeycomb fabric, it facilitates air circulation to prevent possible irritation.
- On its inner surface, it has a hook closure to fasten it to the velour of the cushion.
- Built-in thumb loops support your hand to prevent fatigue.
- Supports the arm and carries the weight across the back and shoulder.
- Universal design fits right or left arm.
- Adjustable straps provide a customisable fit

- Comfortable lightweight material.

### Table of Size:

Measure the length of forearm.



<b>S</b>	20.3-25.4cm (8-10")
<b>M</b>	25.4-30.5cm (10-12")
<b>L</b>	30.5-35.6cm (12-14")
<b>F</b>	20-35.6cm (8-15.5")

### Materials:

TDI Foam, Elastane, PE, Nylon, Polyester.



# Shoulder Support



## Shoulder Support

Warmth & Compression  
Code: 7101

### Indications:

- Shoulder instability
- Rotator cuff tendinitis
- Bursitis

### Features:


- Constructed from heat therapeutic neoprene, provides therapeutic warmth, protection and support for the shoulder and upper arm and aids healing.
- Relief of discomfort and instabilities due to chronic shoulder conditions.
- Helps support injured, weak or arthritic shoulders during occupational and sporting activities and minimizes the potential for future injuries.
- Adjustable hook and loop closure ensures proper fit and compression.
- Universal design fits right or left shoulder.

### Materials:

Neoprene, Cotton, Nylon, Polyester.

### Table of Size

Measure the circumference of the upper arm.



S	22-26 cm (8.66-10.24")
M	26-30 cm (10.24-11.81")
L	30-34 cm (11.81-13.39")
XL	34-38 cm (13.39-14.96")
2XL	38-42 cm (14.96-16.54")
3XL	42-46 cm (16.54-18.11")

# Elbow Support



## Elbow Support

Targeted Relief  
Code: 8101



### Indications:

- Lateral epicondylitis (Tennis elbow)
- Medial epicondylitis (Golfer's elbow)
- Olecranon bursitis
- Weak elbow

### Features:


- Permits limited motion to joint following injury or surgery to avoid stiffness.
- Provides full support of elbow region in conjunction with reducing the insertion of the extensor Carpi Radialis tendon.
- For prevention and control of bursitis and tendinitis, and relief of pain associated with arthritis.
- Adjustable compressive strap reduces tension in wrist and finger extensors during gripping.

### Materials:

Neoprene, Cotton, Nylon, Polyester, PE.

### Table of Size:

Measure circumference of forearm 5 cm (2 inches) below elbow joint.



S

21-23cm (8.27-9.06")

M

23-25cm (9.06-9.84")

L

25-27cm (9.84-10.63")

XL

27-29cm (10.63-11.42")

2XL

29-31cm (11.42-12.20")

3XL

31-33cm (12.20-12.99")

# Elbow Support



## Tennis/Golf Elbow Strap

Targeted Relief  
Code: 8102



### Indications:

- Lateral epicondylitis (Tennis elbow)
- Medial epicondylitis (Glofer's elbow)

### Features:


- The pressure point offers extra pressure on the muscles in the forearm.
- Ideal for any activity involving a strong grip or strain on the forearm and elbow
- Designed to help provide relief from pain and soreness in the forearm and elbow.
- Adjustable hook and loop strap for greater tension if needed.
- Fits left or right.
- One size fits all

### Materials:

Neoprene, Cotton, Nylon, Polyester, PE.

### Table of Size:

Measure circumference of forearm 5 cm (2 inches) below elbow joint.



F

21-34.5 cm(8.27-12.99")

# Wrist Support



## Wrist Band

Adjustable  
Code: 9102



### Indications:

- Wrist sprain, strain.
- Wrist instability.
- Weak wrist.

### Features:


- The pressure point offers specific extra pressure on the wrist.
- Adjustable hook and loop strap insures comfortable fit and stabilization during activities.
- Provides firm support around weak or strained area.
- Heat therapeutic neoprene helps warm muscles and joints.
- The product is a Universal Size (one size fits most), provides adjustable support, fits the Left or Right wrist.

### Materials:

Neoprene, Cotton, Nylon, Polyester, PE.

### Table of Size:

Measure circumference of your wrist.



<b>S</b>	12-14cm (4.72-5.51")
<b>M</b>	14-16cm (5.51-6.29")
<b>L</b>	16-18cm (6.29-7.08")
<b>XL</b>	18-20cm (7.08-7.87")
<b>2XL</b>	20-22cm (7.87-8.66")
<b>F</b>	12-20cm (4.72-7.78")



# Thumb Support



## XR Thumb Brace

Lateral  
Thumb Stay  
Code: 9104



### Indications:

- De Quervain's disease
- Thumb sprain/ strain
- CMC joint injury
- Thumb pain
- Arthritis

### Features:

- Malleable stay holds thumb in desired position for individualized comfort and improved stability.
- Adjustable hook and loop closure ensures comfortable stabilization during activity.
- Additional strap for stabilization during activities.
- Relieves pain and stabilizes CMC joint.
- Supports the wrist and thumb, leaving the palm free.
- Ideal for sprains of thumb, ligament support, and minor fractures of the thumb.
- Neoprene material for maximum comfort in a universal design for left or right.
- The stay can be removed as healing progresses.


### Materials:

Neoprene, Cotton, Nylon, Polyester, PE, Aluminium.

Cured | [www.curedmedical.com](http://www.curedmedical.com)

### Table of Size:

Measure circumference of your wrist.



<b>S</b>	12-14cm (4.72-5.51")
<b>M</b>	14-16cm (5.51-6.29")
<b>L</b>	16-18cm (6.29-7.08")
<b>XL</b>	18-20cm (7.08-7.87")
<b>2XL</b>	20-22cm (7.87-8.66")
<b>3XL</b>	22-24cm (8.66-9.45")
<b>4XL</b>	24-26cm (9.45-10.24")
<b>F</b>	12-22cm (4.72-8.66")


# Wrist Support



## Wrist Brace

**Code: 9108**  
Stabilized Wrist Brace  
(w/Palmar Stay)

**Code: 9109**  
Extreme Stabilized  
Wrist Brace  
(w/Dorsal & Palmar Stays)



### Indications:

- Carpal Tunnel Syndrome
- Colles' fracture
- Early cast removal
- Post-traumatic use

### Features:

- Palmar stay in the 9108 model and palmar and dorsal stays in the 9109 model with cushion pad limit wrist and hand motion for partial immobilization while allowing maximum hand function and provide support weak or injured wrist.
- Contoured aluminium palmar stay is malleable and can be positioned as desired.
- The aluminium stay can be removed as healing progresses.
- Stay can be relocated in "Both hand" designs to fit either left or right hand.
- The adjustable strap provides additional compression around the wrist or hand and ensure proper fit.
- The heat therapeutic neoprene helps warm muscles and joints.
- The compact and effective design allows for easy use with ultimate support and protection. This product has a soft and comfortable construction.


- Free Size version available.

### Materials:

Neoprene, Nylon, Polyester, PE, Aluminium


### Table of Size:

Measure circumference of your wrist.



<b>F</b>	12-22cm (4.72-8.66")
<b>S</b>	12-14cm (4.72-5.51")
<b>M</b>	14-16cm (5.51-6.29")
<b>L</b>	16-18cm (6.29-7.08")
<b>XL</b>	18-20cm (7.08-7.87")
<b>2XL</b>	20-22cm (7.87-8.66")
<b>3XL</b>	22-24cm (8.66-9.45")
<b>4XL</b>	24-26cm (9.45-10.24")

# Wrist Support



## Wrist Brace w/Thumb Spica

**Code: 9110**  
Stabilized Wrist Brace w/Thumb Spica (w/Palmar & Thumb Stays)

**Code: 9111**  
Extreme Stabilized Wrist Brace w/Thumb Spica (w/Dorsal, Palmar & Thumb Stays)

### Indications:

- Carpal tunnel syndrome
- De Quervain's syndrome
- Colles' fracture
- Early cast removal
- Post-traumatic use

### Features:

- Palmar and thumb stays in the 9110 model and palmar, thumb and dorsal stays in the 9111 model with cushion pad limit wrist, thumb and hand motion for partial immobilization while allowing maximum hand function and provide support for weak or injured wrist and thumb.
- Thumb stay for additional support and immobilization.
- Contoured aluminium palmar and thumb stays are malleable and can be positioned as desired.
- The aluminium stays can be removed as healing progresses.
- Stays can be relocated in the "Both hand" model to fit either the left or right hand.
- The adjustable strap provides additional compression around the wrist and/or hand and ensures proper fit.
- The heat therapeutic neoprene helps warm

muscles and joints.


- The compact and effective design allows for easy use with ultimate support and protection.
- Soft and comfortable construction.
- Free Size version available.

### Materials:

Neoprene, Nylon, Polyester, PE, Aluminium

### Table of Size:

Measure circumference of your wrist.



<b>F</b>	12-22cm (4.72-8.66")
<b>S</b>	12-14cm (4.72-5.51")
<b>M</b>	14-16cm (5.51-6.29")
<b>L</b>	16-18cm (6.29-7.08")
<b>XL</b>	18-20cm (7.08-7.87")
<b>2XL</b>	20-22cm (7.87-8.66")
<b>3XL</b>	22-24cm (8.66-9.45")
<b>4XL</b>	24-26cm (9.45-10.24")

# Finger Support



## Finger Splint (Short)

Used in Hospitals  
Code: 10101



### Indications

- Mallet finger
- After treatment
- Postoperative fixation
- Distortion in finger joints

### Features:


- Contoured splints are malleable.
- The splint padded with soft material.
- The splint is easy to wear and remove.

### Materials:

Neoprene, Nylon, Polyester, Metal.

### Table of Size

Measure by laying your hand flat on the table and measuring the height of the finger; the length is from the first knuckle (closest to hand) to the fingertip.




<b>XS</b>	<b>4cm (1.57")</b>
<b>S</b>	<b>5cm (1.96")</b>
<b>M</b>	<b>6cm (2.36")</b>
<b>L</b>	<b>7cm (2.75")</b>
<b>XL</b>	<b>8cm (3.14")</b>

# Lower Extremity

The lower extremity refers to the part of the body from the hip to the toes. The lower extremity includes the hip, knee, and ankle joints, and the bones of the thigh, leg, and foot.

Many people refer to the lower extremity as the leg. In fact, the leg is the part of the body between the knee and ankle joints. The proper way to describe the lower limb is the lower extremity.



# Pelvic Orthosis



## T-POD 6X Pro Trauma Pelvic Orthotic Device Pro

Pulley System 6X  
Used in EMS &  
hospital.  
General: 11301  
Economic: 1130101

1:6 Pulley  
System

6X



General-  
11301



Economic-  
1130101



### Indications

- Pelvic injury
- Pelvic stabilization
- Suspected pelvic fractures
- Internal pelvic hemorrhage
- Pelvic fractures and dislocations

### Features

- T-POD Pro reduced the separation of the normally joined pubic bones by over 70% compared to bed sheets at 44%.
- Provides a quick, safe, and effective method for the initial treatment of pelvic injury and possible pelvic fractures.
- Provides a quick, safe, and effective method for the initial treatment of pelvic injury and possible pelvic fractures.
- Prevent additional damage to veins, arteries and nerves in the area of the fracture.
- Prevent additional movement of the bones near the fracture site.
- Decreases need for administration of blood or blood by-products.
- Reduces patient pain levels and need for pain medication.
- Decreases blood loss or hemorrhage.
- Increases critical mean arterial pressure.
- Reduces the length of hospital stays.


- Remain in place during all forms of patient transport.
- Lowers rate of mortality.
- One person can apply in the field.
- The guide strap controls the circumferential compression by helping to calibrate the gap.
- Small & Lightweight Design
- One Size Fits All
- 100% Radiolucent
- Modulated Compression
- Over-tightening Prevention

### Table Of Size

One size fits all pelvis in the 18.5-56.29" (50-143 cm) hip circumference. Measure around the hip.

### Materials:

Cotton, Nylon, Polyester, Polyethylene




Compression  
Adjustment  
Guide

Application  
History  
Label

Mechanical  
Advantage  
Pulley System

Pull Tab

# Knee Support



## Patella Tendon Strap

Adjustable Fit  
Code: 13101



### Indications:

- Chondropathy and other painful affections of femoropatellar joint.
- Slight patellar instabilities.
- Jumper's knee.
- Mild instabilities of the knee joint.

### Features:


- Wrap around style is easy to put on and take off.
- Trimmed edges and smooth seams.
- Soft Neoprene Blend.
- Use on right or left knee.
- Extra pad provides compression on target area.
- Helps relieve knee pain, supports and stabilizes the kneecap.

### Materials:

Neoprene, Cotton, Nylon, Polyester, Metal

### Table of Size:

The CURED Patella Tendon Strap is a one size fits most knee support and will fit the knee circumference from (30-63 cm). Measure around the knee joint.



# Knee Support



## ADJ Knee Support

Adjustable Fit  
Code: 13104



### Indications:

- Weak knee
- Knee strain/ sprain
- Patellar instability
- Arthritis

### Features:


- Open patella to relieve pressure on the kneecap.
- Neoprene kneecap buttress protects and stabilises the Patella.
- Wrap-around design
- Trimmed edges and smooth seams
- Soft Neoprene Blend
- Easy to put on and take off
- Use on right or left knee
- in two version: Short - Long

### Materials:

Neoprene, Cotton, Nylon, Polyester.

### Table of Size:

Measure around the mid-section of the knee, while in a standing position.



**S** 24.75-31.6(9.74-12.44")

**M** 28.6-38.6(11.25-15.20")

**L** 33.55-44.6(13.2-17.56")

**XL** 38.5-50.6(15.15-19.92")

**2XL** 43.45-56.6(17.1-22.28")

**3XL** 48.4-62.6(19.05-24.65")

**F** 28.6-43.45(11.25-17.1")

# Knee Support



## Knee Stabilizer

Medial & Lateral  
Spiral Stays  
Code: 13106

### Indications:

- Arthritis
- Chondromalacia
- Patellar instability
- Knee strain/sprain
- Post-operative

### Features:


- Open patella to relieve pressure on the kneecap
- Neoprene kneecap buttress protects and stabilises the Patella.
- Adjustable straps offer controlled compression.
- Two medial and lateral spiral stays on both sides of knee providing firm support.
- Retains body heat to warm knee joint and increases blood circulation.
- Trimmed edges and smooth seams
- Use on right or left knee.

### Materials:

Neoprene, Cotton, Nylon, Polyester, Metal, PE

### Table of Size:

Measure around the mid-section of the knee, while in a standing position.



S	29-32 cm (11.41-12.59")
M	32-35 cm (12.59-13.77")
L	35-38 cm (13.77-14.96")
XL	38-41 cm (14.96-16.14")
2XL	41-44 cm (16.14-17.32")
3XL	44-47 cm (17.32-18.50")
4XL	47-50 cm (18.50-19.68")
5XL	50-53 cm (19.68-20.86")

# Knee Support



## ADJ Knee Support (Free Size)

Adjustable Fit  
Code: 131061



### Indications:

- Arthritis
- Chondromalacia
- Patellar instability
- Knee strain/sprain
- Post-operative

### Features:


- Open patella to relieve pressure on the kneecap
- Neoprene kneecap buttress protects and stabilises the Patella.
- Adjustable straps offer controlled compression.
- Two medial and lateral spiral stays on both sides of knee providing firm support.
- Retains body heat to warm knee joint and increases blood circulation.
- Trimmed edges and smooth seams
- Use on right or left knee.

### Materials:

Neoprene, Cotton, Nylon, Polyester, Metal, PE

### Table of Size:

This product is two size that fits most knee support and will fit.



**M**

29-44cm (11.4-17.3")

**L**

44-53cm (17.3-20.86")

# Knee Support



Medial & Lateral Spiral Stays



## Knee Stabilizer (Econ.)

Medial & Lateral Spiral Stays  
Code: 131062

### Indications:

- Arthritis
- Chondromalacia
- Patellar instability
- Knee strain/sprain
- Post-operative

### Features:


- Open patella to relieve pressure on the kneecap
- Neoprene kneecap buttress protects and stabilises the Patella.
- Adjustable straps offer controlled compression.
- Two medial and lateral spiral stays on both sides of knee providing firm support.
- Retains body heat to warm knee joint and increases blood circulation.
- Trimmed edges and smooth seams
- Use on right or left knee.

### Materials:

Neoprene, Cotton, Nylon, Polyester, Metal, PE

### Table of Size:

Measure around the mid-section of the knee, while in a standing position.



S	29-32 cm (11.41-12.59")
M	32-35 cm (12.59-13.77")
L	35-38 cm (13.77-14.96")
XL	38-41 cm (14.96-16.14")
2XL	41-44 cm (16.14-17.32")
3XL	44-47 cm (17.32-18.50")
4XL	47-50 cm (18.50-19.68")
5XL	50-53 cm (19.68-20.86")

# Knee Support



## Knee Immobilizer

Adjustable Fit  
Used in Hospitals  
General: 13301  
Economic: 1330102



### Indications

- Pre-operative
- Post-operative
- Post-traumatic
- Muscle or ligament injuries
- Condition after removal of plaster cast
- Acute treatment after conservatively healed injuries of collateral ligaments
- Undisplaced fractures in the area of knee joint.

### Features

- Knee brace for immobilisation of the knee joint.
- 4-piece with patella Pad.
- Median and lateral splints, which can be freely positioned.
- For children too.

### Materials:

Nylon, Polyester, Metal.


### Table of Size:

The CURED Knee Immobilizer is One Size that Fits Most knee support and will fit the knee circumference from 30 cm to 63 cm. Measure around the knee joint.

Available in two length


Short length 38 cm = Medium Size (for kids)

Long length 50 cm = Large Size (for adult)





# Shin Support



## Hydrotherapy Shin Support (Water\_Resistance Sock)

Hydrotherapy  
Code: 12101



### Indications

- Strengthen leg muscles (quadriceps and hamstrings) and ankle.
- Accelerate the recovery process of athletes after sports injuries in the foot area.
- Effectiveness in water exercise and building endurance in the leg muscles.
- Significant reduction during the treatment.
- There is a safer way to lose weight for obese people by burning more calories, while it may reduce your chances of doing exercise.
- Helps speed recovery in humans with knee pain and rheumatic pain.


### Features:

- Lightweight and easy to use.
- Creating uniform pressure and reliable control over the foot limb, and not creating uncontrollable pressure points during use.
- Leg massage therapy with water flow between the product, wings and calf muscles.
- Significant reduction in the risk of injury compared to similar leg muscle strengthening products such as Weight Calf for injured elderly and overweight people.

- Moving in a water environment because of low gravity is easier, meanwhile it strengthens the muscles.
- The buoyancy of the water helps people balance more. Repetition of movement exercises in water with this product increases strength and endurance and improves flexibility. Besides, exercising in the water has been easier than on land.

### Table of Size

Measure the circumference of the shin.



S	31-34 cm (12.2-13.39")
M	34-37 cm (13.39-14.57")
L	36.5-40 cm (14.37-15.75")
XL	39.5-43 cm (15.55-16.93")
XXL	42-46 cm (16.54-18.11")

### Materials:

Neoprene, PE, Elastane, Nylon, Polyester.

# Ankle Support



## Ankle Support w/Strap

Extra  
Compression  
Code: 14103



### Indications:

- Ankle-joint instabilities
- Ankle strains/ sprains
- Weak or arthritic ankles
- Recurrent ankle distortions

### Features:


- Prophylactic appliance in cases of high-risk of ankle-joint injury.
- Open heel design allows for a comfortable fit and applies less pressure on the sensitive Achilles tendon area.
- Detachable strap may be removed when additional compression is not needed.
- Flexible spiral springs act like tendons to reinforce the ankle.
- Provides comfortable compression for weak or overstressed ankles.
- Relieves pain due to strains or stiffness.
- Free size version available.

### Materials:

Neoprene, Cotton, Elastane, Nylon, Polyester.

### Table of Size:

The Measure circumference just above ankle bone.



F 16-28 cm (6.29-11.02")

XS 16-18 cm (6.29-7.08")

S 18-20 cm (7.08-7.87")

M 20-22 cm (7.87-8.66")

L 22-24 cm (8.66-9.45")

XL 24-26 cm (9.45-10.24")

2XL 26-28 cm (10.24-11.02")

# Foot



15110 (20 mm)



151102 (12 mm)

## Open Toe Cast Sandal

General: 15110  
Economic: 151102

### Indications


- Cast protection
- For early weight bearing and walking on casting.

### Features:

- Features open toe and heel for comfort and breathability.
- Made of ribbed sole to prevent slipping during use.
- Kids versions available.

### Materials:

Nylon, Polyester, Foam, PP.



### Table of Size:

Based on your shoe size, choose the right size according to the table below.

### ADULTS SIZING CHART

<b>XS</b> 29-32 (EU)	<b>L</b> 38-41 (EU)
<b>S</b> 32-35 (EU)	<b>XL</b> 41-44 (EU)
<b>M</b> 35-38 (EU)	<b>2XL</b> 44-47 (EU)

# Sternum Support



## Sternum Support

Used in  
Hospitals  
Code: 2201



### Indications

- Rib injuries
- Chest pain reliever
- Intercostal muscle strain
- Open-heart surgery after sternotomy

### Features:


- Provides compression and support for chest and rib injuries.
- To limit painful chest expansion in case of rib fracture.
- Specifically indicated for open-heart surgery after sternotomy.
- Ideal for support and minimizing movement following rib and sternum fractures.
- Supports soft tissues contusions of the upper chest, sternocostal joint inflammation.
- Pain on breathing or coughing is reduced thereby facilitating lung and breathing physiotherapy.
- It features a full-width velcro closure, adjustable for girth closure with shoulder strap, to prevent slipping.

### Materials:


Cotton, Nylon, Polyester, Polyurethane.  
P50

### Table of Size:

Measure around the chest at the base of the sternum.




A = Men: Circumference of chest measured above breast



B = Women: Circumference of chest measured below breast

<b>XS</b>	61-71 cm (24-28")	<b>XL</b>	102-117 cm (40-46")
<b>S</b>	71-81 cm (28-32")	<b>2XL</b>	117-132 cm (46-52")
<b>M</b>	81-91 cm (32-36")	<b>3XL</b>	132-147 cm (52-58")
<b>L</b>	91-102 cm (36-40")	<b>F</b>	71-147 cm (28-58")



# Bandage



## Elastic Bandage

Medium  
Compression

- Code: 17201 (5cm)
- Code: 17202 (7.5cm)
- Code: 17203 (10cm)
- Code: 17204 (15cm)
- Code: 17205 (20cm)



### Indications:

- General surgical, orthopedic and sports injuries.
- Extremely convenient as a pressure dressing and for skin grafts.
- Can be used for sprains, aches dislocations, painful joints, veins, cramps, skin injuries to tendons and muscles
- Very useful as a medium compression bandage for muscular support


### Sizing:

These bandages are offered in widths of 5, 7.5, 10, 15 and 20 cm, and each bandage has 4 Yards stretched length, equivalent to 3.65 meters.

### Materials:

Polyester, Elastane.

# N95/99 Respirator



## N95/99 Particulate Respirator

W/Valve

N95 Code 5L: 18104

N95 Code 6L: 18106

N99 Code: 18202



### Indications

- Filtration Efficiency >95%/99%
- Standard layers designed
- Lightweight and comfortable
- Epidemic prevention barrier
- Effectively blocks dust, air pollutant and droplet

### Features:

- Can filter out 95/99% of very small (0.3 microns) particles.
- Adaptable nose bar, soft and comfortable to wear.
- Ideal for people at risk of respiratory diseases.
- Designed to fit tightly, creating a seal between your face and the respirator.
- Flat fold design optimizes storage and portability.
- Individual packaging helps prevent contamination.
- Highest level of fluid resistance.
- Helps reduce wearer exposure to certain airborne particles.
- This respirator can be used to protect you from harmful dust particles, PM2.5, vehicle exhaust, woodwork, pollen allergies, and fumes.

- On the inhale airtightly, it cut off pollution into the air.
- Exhalation valve will be opened easily during exhaling and will reduce expiratory air resistance and heat and moisture.
- The aerodynamic design of the air valve guides airflow away from glasses.
- Two-strap design with welded dual point attachment helps provide a secure seal
- Cushioning nose foam Lightweight construction promotes greater worker acceptance and may help increase wear time.


### Table of size:

Available in three sizes (Small, Medium, and Large sizes) depending on face size.

### Materials:

PP, PE, Elastane, Nylon, Polyester.

# N95 Respirator



FILTRATION  
EFFICIENCY >95%



31 PCS N95  
Respirator Filters



Rechargeable  
Fan



## N95 Electronic Respirator

W/Rechargeable Fan  
Code: 18301(General)  
Code: 18302(Sport)



### Indications

- Filtration efficiency >95% for N95.
- Shortness of breathing problems.
- Healthy and comfortable breathing experience.
- Epidemic prevention barrier.

### Features:

- Smart N95 Respirator w/Rechargeable Fan contains miniaturized fresh air module. Deeply purified fresh air is injected into the respirator by the miniaturized fresh air module with proprietary technology. The miniaturized fresh-air module is a comprehensive solution to the traditional respirators with shortness of breath problems, etc. It brings you a healthy and comfortable breathing experience.
- Available Cover in two designs, General and Sport.
- Removable N95 Respirator Filters
- 30 PCS N95 Respirator Filters
- Rechargeable Fan
- Washable Cover

### Table Of Size

One size fits all faces in General and Sports Version.

Cured | [www.curedmedical.com](http://www.curedmedical.com)

### Technical Information:

Rated voltage : DC5V  
Rated power : 0.4w  
Battery capacity : 500mAh  
Battery duration : 5-10 hours  
Weight: 50.5g  
Module size : 75\*37\* 12.5mm  
Executive standard : IEC 60601-1  
:2005+A1:2012  
Effective date : One year.

### Materials

PP, PE, PC, Elastane, Nylon, Polyester.

### Sport Version:




# Compression Socks



## Anti Embolism Stockings

15 - 20 mmHg  
micro fiber  
code: 16312




### Features:

- 18mmHg compression.
- Help reduce the risk of DVT (Deep Vein Thrombosis)
- Pre and post-operative aid for leg discomfort, tired, aching legs, and swollen ankles and feet.
- Aid circulation to help symptoms of mild edema and varicose veins.
- Graduated Compression: maximum compression at the ankle, gradually decreases up the leg, helping aid circulation.

- Microfiber material for comfort, softness and breathability.
- Comfort top band, reinforced toe and a shaped heel cup for extra comfort, durability and best fit.
- Contains 1 Pair.

### Table Of Size

Measure the circumference of each part of your foot according to its shape Size up if you are between sizes.



	S	M	L	XL	2XL	3XL	4XL	5XL
cG	40-44 cm	44-50 cm	47-55 cm	50-58 cm	53-62 cm	55-65 cm	58-68 cm	63-75 cm
cC	25-29 cm	29-33 cm	31-35 cm	34-38 cm	36-41 cm	38-43 cm	41-45 cm	42-48 cm
cB	14-18 cm	18-20 cm	20-22 cm	22-24 cm	24-26 cm	26-28 cm	28-30 cm	30-32 cm



# Compression Socks



## Medical Compression Socks

15 - 20 mmHg  
**LIGHT**  
 20 - 30 mmHg  
**MODERATE**  
 30 - 40 mmHg  
**FIRM**  
 Code: 16201

### Features:

- Medically designed for ultimate comfort and performance.
- Help restore tired, aching legs and swollen ankles and feet.
- Aid circulation to help symptoms of mild edema and varicose veins and as a pre-and/ or post-surgery aid.
- Graduated Compression: maximum compression at the ankle, gradually decreases up the leg, helping aid circulation.
- Microfiber material for comfort, softness and breathability.
- Comfort top band, reinforced toe and a shaped heel cup for extra comfort, durability and best fit.
- Contains 1 Pair.

### Table Of Size

Measure the circumference of each part of your foot according to its shape. Size up if you are between sizes.

	S	M	L	XL	2XL	3XL	4XL
cG	32 - 57	45 - 60	48 - 63	51 - 67	55 - 70	59 - 73	63 - 79
cF	35 - 45	36 - 50	41 - 56	44 - 59	48 - 64	51 - 67	54 - 70
cE	31 - 37	32 - 40	36 - 43	38 - 46	40 - 49	42 - 51	44 - 54
cD	26 - 31	27 - 34	31 - 37	34 - 40	38 - 43	42 - 46	39 - 50
cC	27 - 33	29 - 36	32 - 39	34 - 42	36 - 45	38 - 48	40 - 51
cB1	21 - 24	22 - 27	25 - 31	27 - 34	29 - 36	31 - 37	33 - 39
cB	17 - 19	19 - 21	21 - 23	23 - 25	25 - 27	27 - 29	29 - 32

# Foot



## Halluxflex

Bunion  
Corrector  
Code: 15402



### Indications:

- Correct the deformity of the big toe.
- Pre and post operative.
- Shortness rehab time.

### Features:

- Relieves discomfort and pressure from bunions.
- Reduces risk of bunion progression.
- Designed for both women and men with Hallux Valgus.
- Usable pre or post surgery.
- Innovative knob for precise angle and pressure control.
- Universal design fits left or right foot.
- Soft, breathable medical-grade materials with multiple adjustments.
- All-day comfort, realigns big toe, fits in shoes.
- Adjustable, universal size (one size fits most), unisex.

### Materials:

Nylon, Polyester, Polyurethane, Polyamide, Elastane.


### Table of Size:

One size fits most.

Free Size



# Silicone



## Silicone Heel Cushion

Code: 15403



This product is skin friendly. If you have allergic or sensitive skin, consult with your podiatrist before use.

### Indications

- Plantar Fasciitis
- Heel Pain
- Heel Spur
- Reduces Peak Loading
- Ensures Better Weight Distribution

### Features:

- Dual-density silicone heel cushions reduce shock, relieve pain.
- Washable and transferable between footwear.
- Hypoallergenic.
- Non-toxic.
- Non-flammable.
- Latex free.
- Sold in pairs.
- QTY: 2 pcs

- Moving in a water environment because of low gravity is easier, meanwhile it strengthens the muscles.
- The buoyancy of the water helps people balance more. Repetition of movement exercises in water with this product increases strength and endurance and improves flexibility. Besides, exercising in the water has been easier than on land.

### Materials:

Neoprene, PE, Elastane, Nylon, Polyester.

### Table of Size

Measure the circumference of the shin.

US MEN	US WOMEN	UK	EURO	CURED Size
3 - 4.5	4.5 - 6	2.5 - 4	35 - 37	S
5 - 7	7 - 8.5	5 - 6.5	38 - 40	M
7.5 - 9	9 - 11	7.5 - 9	41 - 43	L
9.5 - 11.5	11.5 - 13	9.5 - 11.5	44 - 46	XL

# Silicone



## Silicone Toe Separator

Four Splints  
Code: 15401



### Indications:

- Hallux valgus syndrome
- POverlapping toes

### Features:

- Helps realign hallux valgus deformities by preventing the big toe from overlapping onto the second.
- Prevents and reduces pain caused by calluses and ingrown toenails.
- QTY: 2 pcs

### Materials:

Pure Silicone



### Instruction for use:

Place Silicone Toe Separator between the first and the second toe. Remove the product once daily for a few minutes.

### Table of Size:

Based on your foot size, choose the right

### CURED SIZING CHART

<b>S</b>	33 - 36 (EU)
<b>M</b>	37 - 41 (EU)
<b>L</b>	42 - 46 (EU)



# Silicone



## Silicone Full Length Medical

Code: 15405



The Full Length Silicone Shoe Insoles are ideal for anyone who suffers from painful foot conditions such as ball of foot pain and plantar fasciitis. The dual density insoles have softer blue spots in the heel and ball of foot area for additional cushioning for those sensitive areas.

### Materials:

Pure silicone

### Table of Size:

Choose based on shoe size.

CURED SIZE	EURO	UK	US WOMEN	US MEN
S	35 - 37	2.5 - 4	4.5 - 6	3 - 4.5
M	38 - 40	5 - 6.5	7 - 8.5	5 - 7
L	41 - 43	7.5 - 9	9 - 11	7.5 - 9
XL	44 - 46	9.5 - 11.5	11.5 - 13	9.5 - 11.5

### Indications:

- Longitudinal and transverse arch support
- Pressure redistribution
- Shock absorption
- Foot fatigue , prolonged standing , and running

### Features:

- Soft blue spots to absorb excess pressure over the heel and metatarsal head
- Optimal redistribution of pressure.
- Memory effect restoration of the original form after use
- Made of high quality silicone.
- Resistant to sweat and water with simple
- QTY: 2 pcs

# Silicone



## Silicone Toe Separator

Code: 15406



The CURED medical-grade Silicone Heel Care Socks provide moisture to dry, rough, and cracked heels and also help reduce friction between heels and footwear. It's washable, hypoallergenic, non-toxic, non-flammable and latex-free. Sold in pairs.

### Indications:

- Providing moisture to dry, rough, and cracked heels
- Reducing friction between heels and footwear.

### Features:

- Silicone coated on the ankle
- High cushioning coefficient.
- Mineral oil soothes and softens skin
- Reduces pain caused by friction between heel and ankle.
- Long functional life.
- Super soft zones.
- Sold in pairs.
- QTY: 2 pcs

### Materials:

Pure Silicone, Polyester, Polyamide, Elastane.

### Instruction for use:

Wash and dry your feet. Slip it on the ankle, and place the heel cups on the back of your heels. Take them off every few hours for at least 15 min and wash your feet.

### Table of Size:

One Size Fits Most.size.



## CURED SIZING CHART

**F**

Free Size

Elastic



## Elastic Elbow Support

Adjustable Fit  
Code: 8201



### Indications:

- Mild sprains and strains
- Elbow protection against cold
- Elbow joints and bone Support

### Features:

- It is made from quality 2-way stretchable material for long-lasting, comfortable support.
- This elastic elbow support retains heat to relieve pain and helps to avoid the frustration resulting from an injury.
- Reduces pains due to sprain and joint stiffness.
- Supports the elbow against sprain and muscle weakness.
- Accelerates recovery time of injuries due to slipping and falling.
- It relieves internal tissue stiffness and reduces swelling by increasing body temperature and improving blood circulation
- QTY: 1 pcs

### Materials:

Elastane, Nylon, Polyester.

### Table of Size:

Measure circumference of forearm 5cm (2 inches) below elbow joint.



**M**

15-17cm (5.9-6.69")

**L**

17-19cm (6.69-7.48")

**XL**

19-21cm (7.48-8.26")

**F**

17-19cm (6.69-7.48")

# Elastic



## Elastic Wrist Support

Code: 9202



### Indications:

- Mild sprains and strains
- Wrist protection against cold
- Wrist joints and bone Support

### Features:


- It is made from quality 2-way stretchable material for long-lasting, comfortable support.
- This elastic wrist support retains heat to relieve pain and helps to avoid the frustration resulting from an injured or weak wrist.
- QTY: 1 pcs

### Materials:

Neoprene, Cotton, Nylon, Polyester, Metal, PE

### Table of Size:

Measure circumference of your wrist.



<b>M</b>	<b>15-17cm (5.9-6.69")</b>
<b>L</b>	<b>17-19cm (6.69-7.48")</b>
<b>XL</b>	<b>19-21cm (7.48-8.26")</b>
<b>F</b>	<b>17-19cm (6.69-7.48")</b>

# Elastic



## Elastic Wrist & Thumb Support

Adjustable Fit  
Code: 9203



### Indications

- Reduces pain
- Wrist and thumb Support
- Mild sprains and strains

### Features


- It is made from quality 2-way stretchable material for long-lasting, comfortable support.
- This elastic wrist and thumb support retains heat to relieve pain and helps to avoid the frustration resulting from an injured or weak wrist.
- QTY: 1 pcs

### Materials:

Elastane, Nylon, Polyester.

### Table of Size:

Measure circumference of your wrist.



<b>M</b>	15-17cm (5.9-6.69")
<b>L</b>	17-19cm (6.69-7.48")
<b>XL</b>	19-21cm (7.48-8.26")
<b>F</b>	17-19cm (6.69-7.48")

# Elastic



## Elastic Ankle Support

Adjustable Fit  
Code: 14201



### Indications

- Mild sprains and strains
- Wrist protection against cold
- Wrist joints and bone Support

### Features:


- It is made from quality 2-way stretchable material for long-lasting, comfortable support.
- This elasticated ankle support is ideal for weak or unstable ankles, sprains, and twisting injuries and reduces pains due to sprain and joint stiffness by increasing body temperature and improving blood circulation.
- Accelerates recovery time of injuries due to slipping and falling.
- Without limiting the ankle joint range of motion
- QTY: 1 pcs

### Materials:

Elastane, Nylon, Polyester.

### Table of Size

The Measure circumference just above ankle bone.



**M**

**15-17cm (5.9-6.69")**

**L**

**17-19cm (6.69-7.48")**

**XL**

**19-21cm (7.48-8.26")**

**F**

**17-19cm (6.69-7.48")**



# Eye Cover



## Eye Cover

Code: 20301



### Indications:

- After eye surgery
- Relieve migraine headaches
- Decreasing light during rest

### Table of Size:

The Measure circumference just above ankle bone.



Adults



Kids


### Features:

- Very lightweight and comfort
- No pressure on the eyes
- Elastic strap for holding the eye cover on your face

### Materials:

Elastane, Nylon, Polyester.

# Elastic



## Elastic Knee Support

(Closed Patella)

Code: 13201



### Indications:

- Mild sprains and strains
- Knee protection against cold
- Knee joints and bone Support

### Features:


- It is made from quality 2-way stretchable material for long-lasting, comfortable support.
- This elasticated knee support is ideal for weak or unstable knee, sprains, and twisting injuries and reduces pains due to sprain and joint stiffness by increasing body temperature and improving blood circulation.
- Accelerates recovery time of injuries due to slipping and falling.
- Without limiting the knee joint range of motion.
- QTY: 1 pcs

### Materials:

Elastane, Nylon, Polyester.

### Table of Size:

Measure around the mid-section of the knee, while in a standing position.



S	33-39 cm (12.9-15.3")
M	36-40 cm (14.1-15.7")
L	40-48 cm (15.7-18.8")
XL	Over 48 cm (18.8"≤)
F	40-48 cm (15.7-18.8")



Comfortable Protection



Covered with two layers of fabric



## Soft Collar

Cervical Collar

Comfortable Protection & Covered with two layers of fabric  
Code: 12013



### Indications:

- Cervical spine contusion
- Acute wry neck
- Cervical spondylosis
- Cervical strain/sprain
- HIVD (Herniated Intervertebral Disc)
- Post-operative

### Features:


- Cervical collar, anatomically shaped.
- Two layers fabric covered.
- Supports the cervical spine in a neutral position.
- Hook and loop closure allows for easy application.
- Washable with convenient hook and loop closure.
- Available in three heights: 4.5cm (1.7 inch), 5.5cm (2.16 inch), 6.5cm (2.5 inch)

### Materials:

TDI Foam, Elastane, Nylon, Polyester

### Table of Size:

- To measure the height of the neck, eyes should be looking straight ahead then measure the distance between the edge of sternum and lower jaw with the neck in the extension position.
- To measure the circumference of the neck, measure from the circumference of the neck at "Adam's" apple.




**S** 24 - 28 cm (9.4-11.02")

**M** 28 - 32 cm (11.02-12.59")

**L** 32 - 36 cm (12.59-14.17")

**XL** 36 - 40 cm (14.17-15.7")

# Spine Thoracic



## Posture Brace

W/Steel Stay

Two  
Flexible Stays  
Code: 31013



### Indications:

- Early kyphosis
- Poor posture
- Weak scapula muscles
- Round shoulder

### Features:


- Additional flexible stays help reinforce back alignment and proper posture and prevent slouching.
- It can also help ease muscular aches and spasms by supporting the soft tissues and evenly distributing forces over the lumbar area.
- Triangle configuration correctly distributes pressure to the tip of the shoulder.
- Crisscrossing straps provide more support and make the brace easy to adjust also provide a custom fit.
- Slimline, so it can be worn under everyday clothes.

### Materials:

Neoprene, Cotton, Nylon, Polyester, Steel.

### Table of Size:

Circumference of chest measured below the breast.



M

54 - 64 cm (21.2-25.1")

L

64 - 74 cm (25.1-29.1")



# Shoulder Support



## Arm Sling

Breathable Fabric  
Code: 73313



The CURED Kids Arm Sling is comfortable and lightweight support. The sling is useful for reducing unwanted or excessive movements at the shoulder joint and clavicle and can be used for right or left shoulder injury. Useful during rehabilitation and post-operative recovery, the support helps prevent overstrain of ligaments, tendons and muscles including the rotator cuff group. It helps relieve pain from bursitis and impingement syndrome and encourages correct alignment of the joint for healing by stabilizing the forearm.

### Indications:

- Brachial nerve injury
- Dislocated shoulder
- Upper arm fracture
- Post-surgical protection

### Features:


- Built-in thumb loops support your hand to prevent fatigue.
- Supports the arm and carries the weight across the back and shoulder.
- Adjustable shoulder strap with comfort pad provides a customisable fit.
- Universal design fits the right or left arm, and its material is comfortable and lightweight.

### Materials:

Cotton, Polyester, Nylon, Polyethylene.

### Table of Size:

Measure the length of forearm.



<b>S</b>	15-20 cm (5.91-7.87")
<b>M</b>	20-25 cm (7.78-9.84")
<b>L</b>	25-30 cm (9.84-11.81")


# Wrist Support



## Wrist Brace w/Thumb Spica

**Code: 91103**  
Stabilized Wrist Brace w/Thumb Spica (w/Palmar & Thumb Stays)

**Code: 91113**  
Extreme Stabilized Wrist Brace w/Thumb Spica (w/Dorsal, Palmar & Thumb Stays)



### Indications:

- Carpal Tunnel Syndrome
- De Quervain's syndrome
- Tendonitis
- Post cast healing
- Scaphoid injuries
- Wrist strains or sprains

### Features:

- Palmar and thumb stays in the 91103 model and palmar, thumb and dorsal stays in the 91113 model with cushion pad limit wrist, thumb and hand motion for partial immobilization while allowing maximum hand function and provide support for weak or injured wrist and thumb.
- Thumb stay provides additional support and immobilization.
- Contoured aluminum palmar and thumb stays are malleable, can be positioned as desired, and removed as healing progresses.
- Stays can be relocated to fit either the left or right hand.
- Soft and comfortable construction
- The heat therapeutic neoprene helps warm muscles and joints.


- The adjustable strap provides additional compression around the wrist and/or hand and ensures proper fit.
- The compact and effective design allows for easy use with ultimate support and protection.

### Materials:

TDI Foam, Elastane, Nylon, Polyester

### Table of Size:

Measure the circumference of your kid's wrist.



<b>S</b>	10 - 12 cm (3.9-4.72")
<b>M</b>	12 - 14 cm (4.72-5.51")
<b>L</b>	14 - 16 cm (5.51-6.3")



# Wrist Support



## Wrist Brace

**Code: 91083**  
Stabilized Wrist Brace  
(w/Palmar Stay)

**Code: 91093**  
Extreme Stabilized  
Wrist Brace  
(w/Dorsal  
& Palmar Stays)

The CURED Kids Extreme/Stabilized Wrist Brace features removable stay to allow variable stabilization and immobilization of the wrist and hand.

### Indications:

- Carpal Tunnel Syndrome
- Colles' fracture
- Early cast removal
- Post-traumatic use

### Features:


- Palmar stay in the 91083 model and palmar and dorsal stays in the 91093 model with cushion pad limits wrist and hand motion for partial immobilization while allowing maximum hand function and provide support for weak or injured wrists.
- Contoured aluminum palmar stay is malleable, can be positioned as desired, and removed as healing progresses.
- Stay can be relocated to fit either the left or right hand.
- The adjustable strap provides additional compression around the wrist and/or hand and ensures proper fit.
- The heat therapeutic neoprene helps warm muscles and joints.
- The compact and effective design allows for easy use with ultimate support and protection.
- Soft and comfortable construction.

### Materials:

Neoprene, Nylon, Polyester, PE, Aluminium

### Table of Size:

Measure the circumference of your kid's wrist.



<b>S</b>	10 - 12 cm (3.9-4.72")
<b>M</b>	12 - 14 cm (4.72-5.51")
<b>L</b>	14 - 16 cm (5.51-6.3")